


Intellego Curing provides a visual indication of applied ultraviolet energy during UV Curing processes. Unlike bulky and complex radiometry equipment, Intellego Curing dosimeters are ultrathin, flexible and also adhesive so they can be used in tight spaces and on dimensional surfaces to provide an accurate dose measurement with ease.

## **Key Features**

- Thin (0.2 mm), flexible indicators for use in narrow spaces
- Intended for use in curing with UV-C light at 254 nm
- Indicator temperature tolerance: 15-70° C | 59-158° F
- Designed with adhesive backing for placement on any object or surface
- 3-year shelf life
- Single use only
- Packaged in resealable mylar pouches
- Private label and customization available
- Polypropylene material is recyclable
- Store indicators in original packaging at 0-40° C | 32-104° F
- Keep away from direct sunlight



UV254-DC-05

## How to use:

- To use, place an indicator on any surface or material that will cured with UV-C devices (254 nm).
- Make sure the yellow indicator area faces the UV light source.
- After exposure, use the color reference to determine if the desired dose of UV energy (mJ/cm²) was achieved.

**SKU** UV254-DC-050

**UV RANGE** 0 - 50 mJ/cm<sup>2</sup>

**DESCRIPTION** 

2.5" x 3.5" dosimeter card

CONTENTS

28 / pouch



Intellego Curing dosimeters are designed to provide a visual indication of applied UV energy. Store at room temperature and away from direct sunlight. **Do not reuse.** Intellego Curing dosimeters are designed and manufactured in Sweden by Intellego Technologies AB.



For more information, please scan the QR code or visit Intellego-technologies.com